

that of the prior lateral-flow immunoassay method. The sensitivity afforded by the present method is so much greater that whereas the previously reported detection limit of lateral-flow immunoassay was 10^5 spores/mL, the estimated detection limit of the present method is 100 spores/mL.

DPA in a 1:1 complex with Ca^{2+} ions is present in high concentration in bacterial spores, and has not been observed in any lifeforms other than bacterial spores. Hence, DPA is an indicator molecule for the presence of bacterial spores. Fortuitously, DPA is also a classic inorganic-chemistry ligand that binds metal ions with high affinity. When bound to Tb^{3+} ions, DPA triggers intense green luminescence of the ions under ultraviolet excitation. The intensity of the luminescence can be correlated with the number density of bacterial spores per milliliter. Moreover, the concentrations of compounds that could potentially give rise to spurious luminescence are typically much smaller than the concentration of DPA, and the strengths with which they bind to Tb^{3+} are of the order of a millionth of that of DPA, so that the desired luminescence signal appears against a dark background.

The figure summarizes the steps of the

present method. A sample suspected of containing bacterial spores is prepared by suspending raw sample material in an aqueous solution that contains Tb^{3+} ions. A volume of $\approx 100\ \mu\text{L}$ of the sample is placed on a test strip — a nitrocellulose membrane on which species-specific antibodies are bound in an area denoted the sample region (area A in the figure). Capillary action moves the spores along the strips. In the sample region, specific binding of membrane-bound antibodies captures and immobilizes the bacterial spores. Next, the strip is exposed to microwave power to release DPA from the spores. The released DPA binds to the Tb^{3+} ions in the solution. Hence, when the strip is exposed to ultraviolet light, the Tb^{3+} ions luminesce green, signaling the presence of the bacterial spores from which the DPA was released.

At the same time, a volume of $\approx 100\ \mu\text{L}$ of a similar solution containing a known concentration of *Bacillus subtilis*, to be used as a standard, is placed on a similarly prepared, parallel membrane denoted the standard strip, which includes an antibody-coated area designated the standard region (area B in the figure). The standard strip is subjected to the same process as is the test

strip. A combination of green luminescence from the region B and a change in color in regions I of both strips indicates that the assay has worked properly. In that case, the ratio of between the intensity of luminescence in region A and that in region B is proportional to the number density of bacterial spores in the sample. The entire assay can be performed in 10 minutes or less.

This work was done by Adrian Ponce of Caltech for NASA's Jet Propulsion Laboratory. Further information is contained in a TSP [see page 1].

In accordance with Public Law 96-517, the contractor has elected to retain title to this invention. Inquiries concerning rights for its commercial use should be addressed to Intellectual Assets Office

JPL

Mail Stop 202-233

4800 Oak Grove Drive

Pasadena, CA 91109

(818) 354-2240

E-mail: ipgroup@jpl.nasa.gov

Refer to NPO-30469, volume and number of this NASA Tech Briefs issue, and the page number.

Live/Dead Bacterial Spore Assay Using DPA-Triggered Tb Luminescence

A relatively simple procedure can be executed in about 20 minutes.

A method of measuring the fraction of bacterial spores in a sample that remain viable exploits DPA-triggered luminescence of Tb^{3+} and is based partly on the same principles as those described in the immediately preceding article. Unlike prior methods for performing such live/dead assays of bacterial spores, this method does not involve counting colonies formed by cultivation (which can take days), or counting of spores under a microscope, and works whether or not bacterial spores are attached to other small particles (i.e., dust), and can be implemented on a time scale of about 20 minutes.

Like the method of the preceding article, this method exploits the facts that (1) DPA is present naturally only in bacterial spores; (2) when bound to Tb^{3+} ions, DPA triggers intense green luminescence of the ions under ultraviolet excitation; and (3) the intensity of the luminescence can be correlated with the concentration of DPA released from

spores and, thus, with the number density of the spores. It has been found that in the case of a sample comprising bacterial spores suspended in a solution, the DPA can be released from the viable spores into the solution by using L-alanine to make them germinate. It has also been found that by autoclaving, microwaving, or sonicating the sample, one can cause all the spores (non-viable as well as viable) to release their DPA into the solution. When the released DPA binds Tb^{3+} ions in the solution and the sample is exposed to ultraviolet light, the solution luminesces, as described in the preceding article.

Therefore, in this method, one divides a sample into two parts. For the first part, germination is used to release the DPA from the viable spores; for the second part, one of the three other techniques is used to release DPA from all the spores. The intensities of the DPA-triggered luminescence of both parts of the sample are measured. Then the frac-

tion of viable spores is calculated as the ratio between the measured luminescence intensities of the first and second parts of the sample.

This work was done by Adrian Ponce of Caltech for NASA's Jet Propulsion Laboratory. Further information is contained in a TSP [see page 1].

In accordance with Public Law 96-517, the contractor has elected to retain title to this invention. Inquiries concerning rights for its commercial use should be addressed to

Intellectual Assets Office

JPL

Mail Stop 202-233

4800 Oak Grove Drive

Pasadena, CA 91109

(818) 354-2240

E-mail: ipgroup@jpl.nasa.gov

Refer to NPO-30444, volume and number of this NASA Tech Briefs issue, and the page number.